


# Disposable spatulas

Size-related designs to fit perfectly with standard 0.2, 0.5, 1.5/2.0, 15 and 50mL microtubes



## Disposable spatulas

Size-related designs to fit perfectly with standard 0.2, 0.5, 1.5/2.0, 15 and 50mL microtubes

- Ideal for dispensing or sample collection of powders or soft solids
- Will not shatter when used with samples frozen in liquid nitrogen
- Hollow core minimises heat transfer — good for hot and cold liquids
- Antistatic, sterile and eco-friendly versions available
- Autoclavable and complies with FDA standards for food contact
- Made of virgin, food-grade polypropylene
- Disposable
- Twin function spatulas, with a spoon at one end and a hollow scoop or pick at the other

Cat. No.	Description	Use with tubes	L x W	Spoon flare	Scoop or pick (L x H)	Pack qty
15356558	Micro antistatic	0.2 and 0.5mL	140 x 3.5mm	n/a	11mm L tapered @ 15° (V-Tip)	300
15326134	Eco-friendly	15 and 1.5mL	140 x 7mm	25 x 14mm	30mm x 3mm	300
15336558	Standard blue	15 and 1.5mL	210 x 7mm	25 x 14mm	30mm x 3mm	300
15346558	Standard opaque	15 and 1.5mL	210 x 7mm	25 x 14mm	30mm x 3mm	300
15376558	Standard opaque sterile	15 and 1.5mL	210 x 7mm	25 x 14mm	30mm x 3mm	100
15366558	Macro green	50mL	310 x 10mm	42 x 23.5mm	30mm L tapered @ 35°	150